

UNITED SOURCES SOUGHT, INC. modernized kit systems

CAGE CODE - 7V5B3



UNITED SOURCES SOUGHT, INC. MODERNIZED KIT SYSTEMS

★ Made in the USA ★ Manufactured and Assembled by VETERANS The state of the art MMPK materials and High-tech construction keep contents such as antibacterial towelettes moist, powder, cell phones, medication and paper documents dry under a broad spectrum of extreme environmental conditions.

The **MMKP** current configuration comes pre

-loaded and ready to use with:

Moist Antibactirial Toweletts for personal and equipment cleanliness Compact Powder Ball with patented formula helps prevent chafing, heat rash, foot fungus

and odor causing bacteria

MILITARY MULTI-PURPOSE KIT SYSTEM = *MMPK*

The *MMPK* System is the first Patented Kit System designed to meet Military, Readiness and Sustainment Specifications.

The Multi-Purpose Waterproof watertight case and cap are designed to protect and keep any essential personal items and valuable accessories secure from external contamination.

 $\mathbf{+}$

Designed to meet Natick Research Criteria the MMPK is able to provide the lightest weight, watertight, puncture resistant multi-fictional case, able to withstand rigorous 7 year testing. The MMPK weighs only 4 ounces.

THE MMPK SYSTEM EXCEEDS MIL SPEC STANDARDS INCLUSIVE OF:

Pressure and temperature tested from Minus 45 °F ro 158 °F to withstand all shipping and inventory storage environments

RF Welded TPU MMPK Case provides water tight Submersion depths up to 300 feet

WWW.USSKITS.COM

239-221-7262

Cage Code 7V5B3

- Industrial Ziplock seal and closure clip, will not crack fade or lose flexibility
- Dust, dirt and sand proof ★ UV and chemical restant ★ Strength <u>Break tested</u>

USE INSTRUCTIONS



1 The Portable MMPK

comes pre-loaded and

ready to use.





2 Turn Grip Cap 1/8 to the left to access Towelettes. **3** Pull up Star Cap to access dry compartment that contains the Powder Ball.

REFILL INSTRUCTIONS



- **1** Open the Exterior MMPK Case Closer Clip located at the end of the MMPK Case.
- **2** Pull apart interior zip closure and remove the empty Towelette Cello Pack.



3 - WIPES Remove protective seal from new Towelette cello pack and slide new pack into the MMPK case. Close both waterproof zip closure and Exterior Portable Case Closer Clip.



4 - POWDER BALL Pull up Waterproof Star Storage Cap to access dry compartment. Remove used Powder Ball. Replace with new Powder Ball and snap close.

COTS Part # - USS1015 © 2017 United Sources Sought, Inc.

Cage Code 7V5B3